

NMIS Enterprise

Features:

- Secures ePHI data in your data center
- Consolidates all of your sites on one server
- Minimizes hardware costs and installation time
- Supports virtual servers and CITRIX environments
- Uses your facility's MS SQL® license (Standard or Enterprise)
- Creates a separate database for each site
- Backups controlled by your database administrator
- Centralizes record review for managers and Radiation Safety Officers
- Share a common Patient Filecard between all or selected sites

RMIS *(Radiopharmacy Management Information System)* Complete Solution to Managing Your Radiopharmacy

Features:

- Online ordering
- Driver dispatch management
- Comprehensive billing module
- Complete pharmacy management
 - Inventory control
 - Prescription processing
- Maintains compliance with ever-changing regulations

PRP Patient Reporting Program

Features:

- Same day report turn around
- Customizable reports
- Improved efficiency
- Electronic signature

Our goal with PRP was to make every cardiologist happy by allowing them to use the expressions that they are comfortable with.

Run your department more efficiently. The PRP provides the user with all tools needed for accurate and timely standardized Nuclar Cardiology Patient Reports that are ASNC, ICANL, and ACR compliant.

INTERFACES

HL7

Time-saving interface which allows the sharing of data between a radiology or hospital information system and the NMIS and Biodose programs.

- Automatic scheduling prior to patient registration
- Interface to any RIS/HIS system
- Completely eliminate transcription errors
- Avoid redundant data entry
- Compatible with NMIS Enterprise
- Send patient dose information back to RIS/HIS

MEDRAD

Eliminate duplicate data entry while reducing your exposure

Take advantage of the reduced radiation exposure and increased dosing efficiency using the Intego PET Infusion System, while maintaining your records in the NMIS program.

- Schedule your patients for the day in NMIS
- Export the schedule to the Intego system. The Intego system will upload the patient ID directly from NMIS
- Dose your patients for the day using the Intego device
- After dosing all of your patients for the day, export the information back to NMIS

Your dosing records will be completed using the information you entered into Intego. Your patients will be dosed, your inventory will be accounted for and your job is done.

DICOM

Time-saving DICOM Worklist interface which allows the sharing of data between a hospital PACS system and the NMIS or Biodose programs.

- Automatic scheduling
- DICOM Worklist interface to any PACS system
- Completely eliminate transcription errors
- Avoid redundant data entry
- Compatible with NMIS Enterprise

Dose Calibrator Interface

The Dose Calibrator interface will allow you to connect certain model dose calibrators to the NMIS or Biodose programs. This interface will capture activity readings in the Dose Calibrator Quality Control, Draw Dose, Generator Elution and Kit Prep modules eliminating typing errors. The Dose Calibrator Interface will also validate that the correct isotope button is pushed. On newer model dose calibrators it will even change the isotope button for you automatically. Please give us a call to see if your dose calibrator is compatible with our Dose Calibrator Interface.



FOLLOW THE RULES. EXCEED THE STANDARD.

BioRx

Worlds First Real-Time Nuclear Pharmacy Management System

Features:

- True real-time processing
- Control all inventory
- Enterprise solution (optional)
- Fast and accurate system
- Easy to understand data entry screens
- QuickBooks compatibility
- Utilizes touch screen and high speed printers
- Bar code technology for enforcing compliance
- Complete billing module
- Accurate reports for monitoring;
 - Health physics statistics
 - Inventory tracking
 - Billing information
- Online ordering

We at ec² Software Solutions are proud to present BioRx, the latest technology in Nuclear Pharmacy Management systems. BioRx is a revolutionary approach to making your pharmacy fast and accurate. Utilizing real-time data entry BioRx was designed to be precise with every dose, every time. Extensive reporting inside BioRx provides each nuclear pharmacy with accurate reports and statistics. Successful Nuclear Pharmacies use BioRx.

BioTRAX QMS

Quality Management System; encompassing all aspects of a PET manufacturing site.

Features:

- Produces electronic batch records
- Document repository
- Centralizes and enforces employee training
- Enforces equipment certification
- Advanced inventory tracking
- Inventory management and re-order reminders
- Ability to monitor facility tasks

BioTrax QMS provides a comprehensive, centralized database for your PET manufacturing site. This electronic information vault includes everything from storing your SOP's and training materials in a document repository to imposing strict stops for critical steps in the process. The extensive inventory management module is used to quarantine, accept, and associate approved raw materials to batch records. We work directly with you to effectuate the roll-out of the software.

In addition the software:

- Has underwent validation testing
- Can be customized to fit any need

BioPointe

Quick, Accurate, and Comprehensive Tracking of Radioactive Deliveries

Features:

- Eliminate delivery errors
- Real time package tracking
- Driver and vehicle dispatch management
- Proof of loading and delivery times
- Runs on Android devices

With BioPointe, you know exactly where all your doses are in the delivery chain.

NMIS & BioDose

Two Unique Products, One Unique Company

Features:

- Patient tracking
- Health physics
- Online ordering
- Logical interface
- Customizable screens
- Predefined and customizable reports
- Data conversion from other software

NMIS and Biodose are designed to meet the needs of every Nuclear Medicine, Nuclear Cardiology and PET facility in today's health care environment. Our software is a complete solution to managing an efficient, compliant and successful department.

Interfaces:

- DICOM
- HL7
- MEDRAD
- Dose Calibrator



400 Apgar Drive
Suite I
Somerset, NJ 08873
800 851 0025 Phone
732 356 8746 Fax

info@ec2software.com
ec2software.com

