

WELCOME TO




Software Solutions



ec² Software Solutions is the leading name in Radiopharmaceutical tracking software and record keeping in the Nuclear Medicine industry. Our software is design to simplify patient scheduling, dosing, inventory control, receiving and waste management. All this combined with our comprehensive health physics suite enables our customer to remain complaint with their regulatory requirements.

ec² product profile NMIS & BIODOSE



FEATURES:

- Dose tracking from arrival through disposition
- Unused doses report to request credits from the pharmacy
- Online ordering
- HL7/DICOM interface option reduces transcription errors
- Complete Health Physics suite to aid in regulatory compliance
- Patients Statistics to document productivity
- Simplified Dosing Screen
- Injection Confirmation Scan

400 Apgar Drive, Suite I
Somerset, NJ 08873

800 851 0025 Phone
732 356 8746 Fax

ec2software.com
info@ec2software.com

SERVICES:

- World class customer service
- Web training available at your convenience
- Onsite training available
- Hours of operations are Monday-Friday 8AM-6PM EST
- Software updates available online

Constancy Test : New Entry

Dose Calibrator Info

Model Manufacturer DC Uniq ID Tech

Serial No. Location Date Time

Comment
(memo)

Unrequired Reason

Test Info

Ref Source	Window	Pred Act	Unit	Reading	%Off
Ce - 57	Ce - 57	110.296	mCi	110.000	-0.268
Ce - 57	F-18	149.470	uCi	149.000	-0.314
Cs - 137	Cs - 137	143.136	mCi	143.000	-0.095
Cs - 137	Ga-67	96.966	uCi	96.500	-0.068
Cs - 137	Tc-99m	121.702	uCi	121.000	-0.577
Cs - 137	Tl-201	316.624	uCi	316.000	-0.197
Cs - 137	Xe-133	19.789	mCi	19.700	-0.450

OK Cancel Automatic

OPTIONS:

- HL7/DICOM interface
- MEDRAD Interface
- Common Filecard
- Dose Calibrator interface
- Dymo Label Printer
- Wrist Band Printer

WORKSTATION HARDWARE REQUIREMENT:

- Intel Core i5 System
- 4 GB RAM
- 250 GB Hard Drive
- Microsoft Windows 7 or Later

[illegible]

ENTERPRISE MANAGEMENT SERVICES:

- Secured ePHI patient data controlled by the DB Administrator
- Reduce exposure to data loss or breach - avoid HIPAA violations and fines
- Unlimited user access for all member sites for convenient reporting
- Minimizes user downtime during upgrades
- Reduce technologist's workload and eliminate transcription errors
- HL-7 eliminates the redundant typing of patient demographics and scheduling information in NMIS or BioDose

From	To	All Rows	Filter	Clear			
Status	Time	Patient ID #	Patient Name	Exam	Flu, Lat	Ordered	Injection Site
	11:11	000-00-0000	PATIENT 01	Cardiac Stress Imaging	No		
13:30	000-00-0000	PATIENT 01	Cardiac Stress Imaging	No			
15:00	000-00-0000	PATIENT 01	Long Vein	No			
16:00	000-00-0000	PATIENT 01	Long Vein	No			

Warning Disables

- QC Complete - ImageID number 1
- Infusion QC - Lactan 14C #1 - Daily
- Canvex QC - Infusion MS Daily QC
- Canvex QC - Infusion MS Weekly QC
- Recovery - Daily

Attention Required

- Disposal - Return to MY PHARMACY
- Area Medication - NOC MED - Daily
- Area Vial Test - NOC MED - Weekly

Prepwork Monitoring Schedule (H)

- Prepwork Monitoring - TECHNICIAN #1 - Daily

Reports

- Report - Daily
- Murphy QC Report - Monthly

Backup

- Backup status - Daily - Done on 05-01-2015 Computer: MANSERENDINO(2) 03:59 PM PM
- Backup status - Weekly - Done on 05-01-2015 Computer: MANSERENDINO(2) 03:58 PM PM

May 2015

Sun	Mon	Tue	Wed	Thu	Fri	Sat
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31				

Today: 5-01-2015
Fri, May 01, 2015

Save Calculator

Backup

Support

Print Labels

Print Reports

Info

Log

Close

Selection Criteria

Review Data

Close Patient

Order Data

Diagnose

Selection History

Software Administrator

Hot Links

Reports

Available Trains

