

DoseSec²ure

DoseSec²ure is a NEW mobile application developed by ec² Software Solutions that integrates with NMIS and BioDose. DoseSec²ure allows you to confirm the right dose is being administered to the right patient. No ePHI data is stored on the mobile device.

HOW IT WORKS ON NMIS & BIODOSE

- Receive and assign a dose to your patient.
- Print DoseSec²ure patient dose label.
- At the point of injection, scan the patient wristband, scan the patient dose label and the Rx label provided by your Radiopharmacy.
- DoseSec²ure will automatically update the NMIS and Biodose dose records.

WHY DO YOU NEED IT?

- Improves patient safety
- Verifies patient doses prior to administration
- Avoid misadministrations

Ask how our HL7/RIS Interface can streamline this process.

